



AirTec 8 Dynamic Mattress Replacement System



Introducing the Cobalt AirTec 8

- The AirTec 8 Dynamic Mattress System is a complete mattress replacement system, designed to provide pressure care, support and comfort in the home or clinical environment.
- Ideal for both Prevention and Resolution, the AirTec 8 is suitable for use by those who are at low-risk all the way through to high-risk, and when used in conjunction with appropriate clinical care the AirTec 8 Dynamic Mattress Replacement System can meet Pressure Care needs up to Stage 4.
- An easy system to set up and use, the Airtec 8 features 20 x 8"/20cm cells, including 4 static head cells. Each cell features a separate lower compartment within the cell to prevent the user bottoming out, meaning no mattress or underlay is required
- A "whisper quiet" digital pump, featuring simple to use adjustments for weight, mode and cycle times, ensure the pump is extremely user friendly.
- The pump not only offers alternating pressure redistribution but it also has a static mode with "constant low-pressure reduction".
- When in alternating mode, the pressure is periodically redistributed by inflating and deflating cells (a 1 in 2 alternation) in order to mimic natural sleep movement patterns which assist in avoiding long-term pressurisation of skin and tissue.
- In constant static mode, the pressure is two-thirds lower when compared to the same pressure level from alternating mode by the means of immersion and envelopment, providing users with a constant low-pressure option.
- Static and alternating Seat Inflation mode is available

cobalthhealth.com.au

Sydney

34/9 Salisbury Road
Castle Hill NSW 2154
02 8853 1100
sales@cobalthhealth.com.au

Melbourne

23/80 Mills Road
Braeside VIC 3195
03 8521 8444
vicsales@cobalthhealth.com.au



Zipper-All-Around

Protects against the ingress of unwanted fluids, and the cover can be removed for laundering



Slip-Resistant Elastic

Convenient straps to secure the mattress in position and avoid movement



CPR Straps

Straps allow for quick deflation of the mattress during emergencies

Features

- Cover is Nylon PU with cotton lining – Biocompatibility, vapour permeable, water resistant, and flame retardant (CA1117).
- The complete all-around zipper, helps reduce the ingress of unwanted fluids and the cover can be completely removed for easy cleaning.
- Ventilated low-air-loss cells help reduce moisture build-up.
- Emergency CPR Tag
- The system has an integrated CPR, Static/Nursing mode, and Pillow Function.
- An anti-slip elasticated strap keeps the mattress in place and prevents movement.
- The easy-to-use connector allows the user to connect the pump to the mattress quickly and easily.
- An air hose Transport Cap seals the air in the mattress to prevent the mattress going down when being transported or transferred, or can keep the mattress inflated in the case of a power outage.
- Audio and Visual Alarm Functions signal Low Pressure or Power Failure faults.
- Control Panel Lock Out function to prevent accidental changes in patient settings
- Durable 125mm TPU Nylon Cell material.
- Flexible pump hooks to suit most bed footboards, with pump bracket for larger beds an optional extra
- Air hose tubes are lined with a spring system to help prevent unwanted kinks or bends and maintain consistent airflow.
- Pressure range 55-120mm
- Air output 7 – 8L/Min

Product Ordering Details & Dimensions

- | | | |
|---------------------|------------------------|-------------------------------------|
| • M22-AMP-S | Single Bed System | 88cm wide x 198cm long x 20cm high |
| • M22-AMP-KS | King Single Bed System | 105cm wide x 198cm long x 20cm high |
| • Pump Dimensions | | 30cm wide x 13.5cm high x 11cm deep |
| • SWL 200KG | | |
| • 12 month warranty | | |

cobalthhealth.com.au

Sydney

34/9 Salisbury Road
Castle Hill NSW 2154
02 8853 1100
sales@cobalthhealth.com.au

Melbourne

23/80 Mills Road
Braeside VIC 3195
03 8521 8444
vicsales@cobalthhealth.com.au