

PB Deep Back

EXCEPTIONAL SUPPORT AND COMFORT



The Invacare® Matrx® PB Deep Back is contoured to provide increased midline support and positioning. The extensive range of sizes provides exceptional comfort and support for a wide range of individuals with moderate to high lateral trunk positioning needs.



PB Deep Back



Foam:

- Molded foam wraps around shell for maximum comfort and protection
- Includes removable positioning sacral/lumbar support pad

Outer cover:

- Moisture resistant and breathable
- Optional privacy flap is available

Back shell:

- Contoured shell provides exceptional midline support and positioning
- Durable, lightweight ABS back shell includes 7" contour depth and built-in plate for headrest attachment

Hardware:

- "Easy Set" hardware installed and adjusted in minutes - adjustments can be made while user is seated in the wheelchair

Accessories

Chest Harness Bracket



Privacy Flap



Split



Heavy-Duty Easy Set



Extended Range










Fixed Space Saver



Mounting Pins



Tech Specs

	 Height	 Width	 Depth adjustment	 Angle adjustment	 Product weight	 Weight Capacity	 Limited Warranty
PB Deep	12", 16", 20"	14", 16", 17", 18", 20"	2" forwards, 2" backwards	40° range**	From 1.3kg (2.9lb)	136kg (300lb)	2 years
PB Deep Heavy Duty	16", 20"	20, 22, 24, 26, 28, 30"	2" forwards, 2" backwards	40° range**	from 2.7kg (6lb)	226kg (500lb)	2 years

** Limitation of depth/angle adjustment if used on chairs narrower than model width



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®